AUTOCLAVE (150-160 L) Top loading

Specifications:

Capacity	: 150-160L
• Model	: Fully Automatic
Chamber type	: Top loading
• Timer	: 0-10hr
Pressure	: 2.0 Kg/ cm ² at 130 ^o C
Pressure gauge	: 0-4 Kg/ cm ²
Temperature accuracy	:±0.5 º C
• Sterilization Temperature	: Ambient 10 º C -130 º C
• Sterilization time range	: 0-9h
Material of chamber	: SS 304 grade/316Ti grade with corrosion resistance
Pressure resistance	: Chamber pressure should resist 60 psi for safety
Material exterior	: Epoxy resin powder coated steel
Display	: LED/LCD to indicate working status of temp. & pressure
Alarms	: Automatic shutdown with beep reminding after
	sterilization and recycle of operation is completed.
 Lid closing/opening 	: Pedal free system for lid closing
 Safety valves 	: Safety valve for 20 psi with an extra safety valve for 25 psi
Safety features	: Over temperature, Low level water limit, over current, auto -release against over pressure
Baskets	: 2 nos. Of baskets (Mesh/Perforated)
Operating Voltage	: 230 ±10VAC, 50 Hz.
operating voltage	. 200 210

Other Important Requirements

- A control lockout to prevent a cycle start when the door is not closed and locked
- Steam pressure lock to prevent operator from opening if pressure exists in the chamber.
- After reaching normal temperature, door should be opened by a single switch mode, preferably
- Water auto fill for steam generation
- Internal memory for documentation of ~500 sterilization cycles
- Three Years comprehensive warranty with 2 years free AMC
- Local Service Support should be within 24 hours
- Additional features/accessories if any that can potentially increase the productivity and safety of the instrument should be quoted as optional items.